

MEDIRAD Duoburn CD/DVD Disc Publisher

Key Features:

Low cost alternative to duplicating films

Intuitive User interface turns novices into experts in minutes

Receives and stores DICOM images onto a custom labeled CD/DVD coming from a DICOM-network

Touch Screen Ready

Windows 10 compatible

AES 256bit Encryption

Burn studies from your web browser

Burn studies automatically as soon as it is received

Build-in Webserver Service

Disc Spanning for large Studies

Automatically remove Studies after publishing

Easy enable/disable the Web Module

Include new Weasis Viewer

On-demand disc creation and labeling directly from modality workstations

Physicians have access to multiple studies on one single CD/DVD.

Automatically generates a DICOMDIR root file and HTML, and can be stored on each produced CD or DVD

Print dynamic information such as Patient ID, Patient Name, Study Description, Institution etc.

Studies can be conveniently viewed from disc on a workstation using the included DICOM viewer or specialized viewers. Works on Windows XP SP 3, Vista, Windows 7, 8, 8.1 and 10 Meets all industry standards including DICOM part 10, IHE PDI (Dicom Structured Report)

Multiple Viewers

You have the ability to automatically include multiple viewers. Windows, OSX or Linux based.

Additional Features:

Dicom Query/Retrieve from PACS

Patient Study Export to Folder or USB through shared network

