

# PPMAMMO QC & Viewer


## PPMAMMO QC & Viewer

PACSPLUS MAMMO QC is a quality control station automatically generates symmetrical 2 on 1 stitched images. Without complicated hanging protocol of other MAMMO viewers, PPMAMMO QC dramatically increase the speed of image display and storage capacity of typically high-volume of breast images as well as enhance ease of use. PPMAMMO makes it possible to work interactively and intelligently through information sharing and offer fast access to patient images.




## Features

- » Auto trimming process of image increases storage capacity
- » Auto stitching of two images with aligning height of both nipples
- » Reduce volume of the image without image loss
- » W/L, Zoom in/out, magnifying glass, etc
- » Import JPG, BMP, TIFF, etc to DICOM
- » Publisher provides patient CD/DVD
  - DICOM Print
  - Support DICOM storage SCP/SCU




Trimming unnecessary space



Matching and Stitching R/L



Report



Zooming