

PPSCREEN

PPSCREEN (PACSPLUS Screen Viewer)

PACSPLUS SCREEN is a software for multi diagnosis. Image search, display, save, report can be done in PPSCREEN as well as database management function. It is very efficient to interface with HIS to manage patient data exam and image. It is a DICOM compatible and easy to interface with other PACS and modality.



Features


- » Multi exam display
- » Multi exam reading
- » Code reading
- » R/L sync (Zoom, Panning, Magnifying glass)
- » Comparison mode with previous exam and various layout
- » Various image processing and editing
- » DICOM image display and Non-DICOM image (JPEG, TIFF, BMP) to DICOM format file
- » DICOM Storage, Query & Retrieve, DICOM DIR
- » Backup to CD/DVD
- » Interface with other DICOM Viewer



Multi exam display



PPSCREEN Viewer



Display layout



Comparison mode