

# PPFOCUS

## PPFOCUS (PACSPLUS Focus Viewer)


Ophthalmology PACS digitizes the ophthalmology practice. Now with the introduction of PACS technology into the ophthalmology practice, the retinal and fundus images can be digitized providing superior images that can be enhanced and manipulated for optimal viewing.

Workstation with all the necessary functions for Ophthalmology




## Features

- » Window width and level adjustment
- » Zoom in and Zoom out capabilities
- » Various measuring functions (length, angle and etc.)
- » Rotation/Reversal of images
- » Diagnostic opinions on images
- » Convenient thumbnails for inquiring images
- » Various image layout
- » Series layouts on images
- » Indicate previous exam lists
- » Easy compare before/after images
- » Stereo View (Display three dimensional image using 3D spectacles)
- » Worklist
- » Panorama View
- » RGB Split
- » C/D Ratio (Ratio of disc and cup of optic nerve's ending)




Stereo View



RGB Split and Stereo View



Full Image Viewer



C/D ratio | Panorama View