

PPW

PPW (PACSPLUS Viewer)


PPW is the multi-modality viewing software and PPW provides viewing and manipulation of radiological data including images, reports, patient status and clinical information.

An easy to use DICOM Viewer that is used in thousands of hospitals, imaging centers and clinics worldwide. PACSPLUS Workstation incorporates many features required for day-to-day work. Window/level presets, synchronized stacking, reference lines, cine function and integrated CD & DVD burning, as well as advanced capabilities whole spine stitching. PACSPLUS also includes Hanging Protocol builder with relevant priors and Key Images.




Features


- » Various image display functions
- » Thumb-nail display mode
- » Display scout-line of CT, MR and PET, etc
- » Hanging protocol of image display (Image layout, series layout)
- » Comparison of current and historical previous images
- » Cobb angle of Scoliosis (LT, RT)
- » CD/DVD burning (Embedded Auto-run Viewer)
- » Non-Linear Look-Up-Table
- » Long-bone stitching



PACSPLUS Publisher (PPP)



Non-Linear Look Up Table (LUT)



Scout View



Cobb angle of Scoliosis (LT, RT)